

## Acecide®-C High-Level Disinfectant

An easier, safer, more efficient reprocessing solution.



**OLYMPUS**  
**Acecide®-C**  
HIGH-LEVEL DISINFECTANT AND STERILANT  
FOR USE IN OLYMPUS OER-PRO ONLY  
**SOLUTION 1 875mL**  
Cassette Bottle for OER-PRO

**ACTIVE INGREDIENTS:**  
Peracetic Acid ..... 4.8%  
Water & Inert Ingredients ..... 95.2%  
Total ..... 100.0%

**KEEP OUT OF REACH OF CHILDREN—DANGER—PRECAUTIONARY STATEMENTS:**  
Hazard if swallowed. Do not get in eyes, on skin or clothing. Wear protective eyewear (goggles, face shield or safety glasses) and rubber gloves when handling. Wash hands thoroughly with soap and water after handling. Do Not Use for disinfecting contact lenses, contact or lens. **WARNING:** See instructions for Use for Recommended Use of ACECIDE-C (See Solutions).  
**Failure to follow instructions for Use may result in patient or user injury. DIRECTIONS:** In preparation of ACECIDE-C (See Solution), follow the instructions for Use. ACECIDE-C (See Solution) is prepared by mixing 1 part of ACECIDE-C Solution 1, 1 part of ACECIDE-C Solution 2 and 18 parts water. The use must adhere to the instructions for Use since any modification will affect the safety and effectiveness of the peracetic acid. **INTENDED USE:** ACECIDE-C Solution 1 without heat diluting and buffering. Heat the ACECIDE-C solution and suitable Automated Endoscope Reprocessor (AER) treated before use. **INTENDED USE:** ACECIDE-C Solution 1 must be combined with ACECIDE-C Solution 2 and water in the correct proportion before use. ACECIDE-C (See Solution) is intended for the high-level disinfection or sterilization of reusable heat-stable critical and semicritical medical devices (such as flexible endoscopes). ACECIDE-C (High-Level Disinfectant and Sterilant) can be used in Automated Endoscope Reprocessors (AER) systems when approved by the manufacturer of the AER. **USER PRECAUTIONS:** It is a violation of Federal law to use this product in a manner inconsistent with its labeling. The user should be advised that the manufacturer and distributor of medical devices are not responsible for the handling of such substances such as liquid chemical germicides. Additional information about ACECIDE-C (High-Level Disinfectant and Sterilant) can be obtained from Best Solutions, Inc., at 1-888-222-5247. **IN CASE OF A MEDICAL EMERGENCY:** In case of contact with eyes, immediately flush with large amounts of water for least 15-20 minutes. Get prompt medical attention. In case of contact with skin, remove contaminated clothing and wash area with water. Get medical attention if irritation persists. In case of ingestion, do not induce vomiting. Drink large amounts of water immediately to dilute. Seek medical attention immediately. In case of inhalation, move to well-ventilated area and seek medical attention if irritation persists. For chemical emergency, spill, leak, fire, exposure or accident, call CHEMTEL, day or night (800) 424-9300, for District of Columbia or outside the U.S.A. call (702) 527-5887.

**NOTE TO PHYSICIAN:** Possible mucous membrane damage may compromise gastric integrity. **STORAGE:** Store upright in shipping carton. Do not expose to direct sunlight. ACECIDE-C Solution 1 should be stored in its original container or transferred into a clean, opaque (UV-protected) HDPE bottle. ACECIDE-C should not be used after expiration date. **(CONTINUED ON BACK LABEL)**  
Manufacturer: Best Solutions, Inc., 90 8881, 1900 Valley, CA 95044  
Distributor: Olympus America, Inc., 3500 Corporate Parkway, Center Valley, PA 18034-2615  
MADE IN USA A1-105-165-000-000

**OLYMPUS**  
**Acecide®-C**  
HIGH-LEVEL DISINFECTANT AND STERILANT  
FOR USE IN OLYMPUS OER-PRO ONLY  
**SOLUTION 2 875mL**  
Cassette Bottle for OER-PRO

**SOLUTION 2 BUFFER CONCENTRATE**  
**ACTIVE INGREDIENTS:**  
Sodium Phosphate Buffer ..... 4.1%  
Water & Inert Ingredients ..... 95.9%  
Total ..... 100.0%

**KEEP OUT OF REACH OF CHILDREN—DANGER—PRECAUTIONARY STATEMENTS:**  
Hazard if swallowed. Do not get in eyes, on skin or clothing. Wear protective eyewear (goggles, face shield or safety glasses) and rubber gloves when handling. Wash hands thoroughly with soap and water after handling. Do Not Use for disinfecting contact lenses, contact or lens. **WARNING:** See instructions for Use for Recommended Use of ACECIDE-C (See Solutions).  
**Failure to follow instructions for Use may result in patient or user injury. DIRECTIONS:** In preparation of ACECIDE-C (See Solution), follow the instructions for Use. ACECIDE-C (See Solution) is prepared by mixing 1 part of ACECIDE-C Solution 1, 1 part of ACECIDE-C Solution 2 and 18 parts water. The use must adhere to the instructions for Use since any modification will affect the safety and effectiveness of the peracetic acid. **INTENDED USE:** ACECIDE-C Solution 1 must be combined with ACECIDE-C Solution 2 and water in the correct proportion before use. ACECIDE-C (See Solution) is intended for the high-level disinfection or sterilization of reusable heat-stable critical and semicritical medical devices (such as flexible endoscopes). ACECIDE-C (High-Level Disinfectant and Sterilant) can be used in Automated Endoscope Reprocessors (AER) systems when approved by the manufacturer of the AER. **USER PRECAUTIONS:** It is a violation of Federal law to use this product in a manner inconsistent with its labeling. The user should be advised that the manufacturer and distributor of medical devices are not responsible for the handling of such substances such as liquid chemical germicides. Additional information about ACECIDE-C (High-Level Disinfectant and Sterilant) can be obtained from Best Solutions, Inc., at 1-888-222-5247. **IN CASE OF A MEDICAL EMERGENCY:** In case of contact with eyes, immediately flush with large amounts of water for least 15-20 minutes. Get prompt medical attention. In case of contact with skin, remove contaminated clothing and wash area with water. Get medical attention if irritation persists. In case of ingestion, do not induce vomiting. Drink large amounts of water immediately to dilute. Seek medical attention immediately. In case of inhalation, move to well-ventilated area and seek medical attention if irritation persists. For chemical emergency, spill, leak, fire, exposure or accident, call CHEMTEL, day or night (800) 424-9300, for District of Columbia or outside the U.S.A. call (702) 527-5887.

**NOTE TO PHYSICIAN:** Possible mucous membrane damage may compromise gastric integrity. **STORAGE:** Store upright in shipping carton. Do not expose to direct sunlight. ACECIDE-C Solution 2 should be stored in its original container or transferred into a clean, opaque (UV-protected) HDPE bottle. ACECIDE-C should not be used after expiration date. **(CONTINUED ON BACK LABEL)**  
Manufacturer: Best Solutions, Inc., 154 Melton Drive, Walnut, NY 61084  
Distributor: Olympus America, Inc., 3500 Corporate Parkway, Center Valley, PA 18034-2615  
MADE IN USA A1-105-167-000-000

# Acecide®-C High-Level Disinfectant

Acecide-C is a room temperature, peracetic acid-based liquid chemical germicide designed for the Olympus OER-Pro—proven for material compatibility by Olympus for Olympus endoscopes.

## Features and Benefits

- 1 Unique, easy-to-handle cassette bottles allow for quick and easy chemical loading while minimizing exposure to the disinfectant—enhancing safety for the user.
- 2 Fast acting PAA delivers high-level disinfection in 7 minutes resulting in a short 26-minute OER-Pro<sup>1</sup> cycle time.
- 3 Time-saving test strip verifies efficacy in just 13 seconds—changing from white to black for clear interpretation of results.
- 4 Convenient room-temperature (20°C) use protects equipment while saving on utility costs.
- 5 Only PAA available that has been validated for material compatibility with Olympus endoscopes by the endoscope manufacturer.

## Product Overview

Acecide-C is an aldehyde-free, reusable liquid chemical germicide designed for exclusive use with the OER-Pro. Acecide-C is delivered in a set. One set of Acecide-C is loaded into the disinfectant drawer of the OER-Pro, when Acecide-C Solutions 1 and 2 are combined automatically with water to deliver an Acecide-C Use-Solution for high-level disinfection with PAA as the biocidal agent. Acecide-C is reusable for up to 5 days or when the Acecide test strip reveals that the concentration is below the Minimum Recommended Concentration, whichever comes first.

## Product Characteristics

- Minimum Recommended Concentration (MRC): 2,000 ppm (0.2%)
- High-Level Disinfection: 7 minutes at 20°C (68°F)
- Reuse Period: Up to 5 days
- No activation required
- Active Ingredient: Peracetic Acid

## Ordering Information

For more information, contact your local sales representative or call 800-848-9024.

### **ACECIDE-C**

Acecide-C High-Level Disinfectant  
(6 sets of Acecide-C)

### **ACECIDE-C TEST**

Acecide-C Test Strips  
(100 count)

**Note:** Please read and follow the Instructions for Use and MSDS prior to handling Acecide. Minimize exposure and contain or exhaust the fumes from all liquid chemical germicides. Always follow the label directions.

<sup>1</sup> The OER-Pro can be used with Acecide-C to achieve high-level disinfection in 7 minutes, or with ALDAHOL®, a ready-to-use disinfectant that achieves high-level disinfection in as low as 10 minutes.



Olympus and OER-Pro are trademarks or registered trademarks of Olympus Corporation, Olympus America Inc., and/or their affiliates. ALDAHOL is a trademark of Healthpoint, Ltd. Acecide is a registered trademark of Saraya Co., Ltd.